Homework #3

1. Stock Solution of Lysozyme: $10 \text{ prl} \times 50 \text{ mg} = 500 \text{ mg} \text{ or } 0.5 \text{ g of}$ 1 prl lysozyme

Dissolve 0.5 g of lysozyme in about 5 ml of TE. Once dissolved, dilute the stock solution to 10 ml with TE in a graduated cylinder.

For your protein purification:
$$10 \text{ mg/s} \times 1 \text{ ml} = 0.2 \text{ ml of } 50 \text{ mg/ml}$$

 $1 50 \text{ mg/s}$ lysozyme

In my protein purification, I would add 0.2 ml or 200 μ l of 50 mg/ml lysozyme.

2. 1x Tris-glycine

Tris:
$$1.0 \cancel{k} \ge 25 \text{ m/M} \ge 1 \text{ Mote} \ge 121.4 \text{ g} = 3.04 \text{ g of Tris}$$

1 $\cancel{k} = 1000 \text{ m/M} \text{ Mote}$

For 5X: 3.04 g of Tris x 5 = 15.2 g of Tris

Glycine: $1.0 \cancel{k} \ge 250 \cancel{mM} \ge 1$ Mole $x \cancel{75.07 \text{ g}} = 18.77 \text{ g of Glycine}$ 1 \cancel{k} 1000 \cancel{mM} Mole

For 5X: 18.77 of Glycine x = 93.85 g of Glycine

SDS: $0.1 \% = \frac{0.1 \text{ g}}{100 \text{ ml}}$

$$\frac{1.0\cancel{L}}{1} \ge \frac{1000 \text{ mf}}{\cancel{L}} \ge \frac{0.1 \text{ g}}{100 \text{ mf}} = 1 \text{ g of SDS}$$

For 5X: 1.0 g of SDS x 5 = 5.0 g of SDS

Dissolve 15.2 g of Tris, 93.85 g of Glycine and 5.0 g of SDS in about 500 ml of water. Once dissolved, dilute the stock solution to 1.0 L or 1000 ml in a graduated cylinder.

3. 1x Tris-glycine from 5x:

= 0.2 L or 200 ml of 5x Tris-glycine buffer

Dilute 200 ml of 5x Tris-glycine buffer to 1.0 L with H₂O in a graduated cylinder to make 1x Tris-glycine buffer. (Also, you could say 200 ml of Tris-glycine plus 800 ml of H₂O.)

4-8 Check your Biochemistry book for 3 and 1 letter codes of the amino acids.

- 4. Glycine, Alanine, Valine, Leucine, Isoleucine, Methionine, Proline, Phenylalanine and Tryptophan
- 5. Same as #4
- 6. Lysine and Arginine
- 7. Serine, Threonine, Asparagine, Glutamine, Tyrosine, Cysteine, Lysine, Arginine, Histidine, Aspartic Acid, Glutamine Acid
- 8. Aspartic and Glutamic Acid